

60th Medical Group (AMC), Travis AFB, CA
INSTITUTIONAL ANIMAL CARE AND USE COMMITTEE (IACUC)
FINAL REPORT SUMMARY

(Please type all information. Use additional pages if necessary.)

PROTOCOL #: FDG20130030A

DATE: 23 April 2014

PROTOCOL TITLE: Pilot Study: Wound healing with CorMatrix® in the rabbit (*Oryctolagus cuniculus*).

PRINCIPAL INVESTIGATOR (PI) / TRAINING COORDINATOR (TC): Lt Col Antoinette Shinn

DEPARTMENT: 60 MDG/SGSE

PHONE #: 423-7267

INITIAL APPROVAL DATE: 3 June 2013

LAST TRIENNIAL REVISION DATE: N/A

FUNDING SOURCE: SG Office

1. RECORD OF ANIMAL USAGE:

Animal Species:	Total # Approved	# Used this FY	Total # Used to Date
<i>Oryctolagus cuniculus</i>	10	5	5

2. PROTOCOL TYPE / CHARACTERISTICS: (Check all applicable terms in **EACH** column)

<input type="checkbox"/> Training: Live Animal	<input type="checkbox"/> Medical Readiness	<input type="checkbox"/> Prolonged Restraint
<input type="checkbox"/> Training: non-Live Animal	<input type="checkbox"/> Health Promotion	<input type="checkbox"/> Multiple Survival Surgery
<input checked="" type="checkbox"/> Research: Survival (chronic)	<input type="checkbox"/> Prevention	<input type="checkbox"/> Behavioral Study
<input type="checkbox"/> Research: non-Survival (acute)	<input type="checkbox"/> Utilization Mgt.	<input type="checkbox"/> Adjuvant Use
<input checked="" type="checkbox"/> Other (Model/Technique Perfecting)	<input type="checkbox"/> Other (Treatment)	<input type="checkbox"/> Biohazard

3. PROTOCOL PAIN CATEGORY (USDA): (Check applicable) ☐ C ☒ D ☐ E

4. PROTOCOL STATUS:

***Request Protocol Closure:**

☐ Inactive, protocol never initiated
☐ Inactive, protocol initiated but has not/will not be completed
☒ Completed, all approved procedures/animal uses have been completed

5. FUNDING STATUS: Funding allocated: \$4850.00 Funds remaining: \$ 0.00

6. PROTOCOL PERSONNEL CHANGES:

Have there been any personnel/staffing changes (PI/CI/AI/TC/Instructor) since the last IACUC approval of protocol, or annual review? ☒ Yes ☐ No

If yes, complete the following sections (Additions/Deletions). For additions, indicate whether or not the IACUC has approved this addition.

ADDITIONS: (Include Name, Protocol function - PI/CI/AI/TC/Instructor, IACUC approval - Yes/No)

Report Documentation Page		Form Approved OMB No. 0704-0188
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.		
1. REPORT DATE 15 APR 2014	2. REPORT TYPE Final	3. DATES COVERED 3 Jun 2013 - 15 Apr 2014
4. TITLE AND SUBTITLE FDG20130030A "Pilot study: Wound healing with CorMatrix in the rabbit (Oryctolagus cuniculus)."		5a. CONTRACT NUMBER
		5b. GRANT NUMBER
		5c. PROGRAM ELEMENT NUMBER
6. AUTHOR(S) Lt Col Antoinette Shinn, J. Kevin Grayson, DVM, Ph.D.		5d. PROJECT NUMBER FDG20130030A
		5e. TASK NUMBER
		5f. WORK UNIT NUMBER
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Clinical Investigation Facility David Grant Medical Center 101 Bodin Circle Travis AFB, CA 94535		8. PERFORMING ORGANIZATION REPORT NUMBER
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES) Clinical Investigation Facility David Grant Medical Center 101 Bodin Circle Travis AFB, CA 94535		10. SPONSOR/MONITOR'S ACRONYM(S)
		11. SPONSOR/MONITOR'S REPORT NUMBER(S)
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release, distribution unlimited		
13. SUPPLEMENTARY NOTES		

14. ABSTRACT

Objective: Technique perfecting for wound healing model in the rabbit **Results:** Using the five rabbits for technique perfecting proved very beneficial for the investigator. Multiple lessons were learned. 1. The CorMatrix covered the ear wounds best when they were made slightly larger than the wound. Originally, the ear wounds and CorMatrix patches were both cut with a 6mm biopsy punch. This made stitching the patch in place on all sides difficult. The patch was flush with the wound on one side and not necessarily on the opposite side. When the wound was made with a 5mm biopsy punch and the patch with the 6mm punch, the patch easily covered the entire wound on all sides when secured in place. 2. Four wounds per ear was optimum. More than four wound made it difficult (if not impossible) to avoid lacerating a large blood vessel or creating a wound proximal to the ear orifice (an area difficult to place a wound dressing). It also made it difficult to maintain consistent spacing between each wound in an area of the ear that could have a dressing applied. 3. The original spacing between the wounds was 30 mm. That is much larger than necessary or desired. The 30mm distance made applying a dressing to cover all four wounds challenging. Making the wounds 10-15mm apart kept the wounds far enough apart to easily suture the treatment patches in place and not cross contaminate the controls. It also made it much easier to apply a dressing to the ears. 4. Each stitch required a minimum of four throws to tack the CorMatrix in place on four quadrants. The PDS suture has memory and easily unties when not adequately secured. In the first animals, only two throws per stitch were performed and the patches did not remain in place. Upon examination of corresponding histology slides there was no evidence of the CorMatrix patches. When four throws were used to secure the patches, the patch was verified by the histology results. 5. The 5-0 PDS worked much better than 5-0 Vicryl for suturing the patches in place. The Vicryl did not move through the tissue and cartilage smoothly (it was very rough). 6. The Tegaderm dressing worked best when only a small piece was (cut to size) used on the ear. However it would only last one day at best, before coming off the ear. Bacitracin ointment applied twice a day worked well to keep the wound beds moist. However, the bacitracin ointment applied to the wounds with a Vaseline/petroleum gauze place over it appeared to be the best ear dressing (for maintain wound moisture) and would remain in place for at least a day (possibly two days) before changing. 7. EMLA cream could also be applied post operatively before the dressing is applied for postoperative analgesia. 8. The wounds healed very quickly and were completely healed before 28 days post wounding.

15. SUBJECT TERMS

US Air Force, Medical Service, Medical Research, Graduate Medical Education

16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT	18. NUMBER OF PAGES	19a. NAME OF RESPONSIBLE PERSON
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified	UU	3	

None

DELETIONS: (Include Name, Protocol function - PI/CI/AI/TC/Instructor, Effective date of deletion)

Dr. J. Kevin Grayson, Co-Investigator (Amendment 1 – approved, November 8, 2013)

7. PROBLEMS / ADVERSE EVENTS: Identify any problems or adverse events that have affected study progress. Itemize adverse events that have led to unanticipated animal illness, distress, injury, or death; and indicate whether or not these events were reported to the IACUC.

None

8. REDUCTION, REFINEMENT, OR REPLACEMENT OF ANIMAL USE:

REPLACEMENT (ALTERNATIVES): Since the last IACUC approval, have alternatives to animal use become available that could be substituted in this protocol without adversely affecting study or training objectives?

A small animal model using rodents (mice/rats) with a splinted skin wound should be evaluated. In full thickness excision, the mechanical structure of the dermis is completely disrupted, and various forces in the surrounding dermis and the wound site can change (usually decrease) the dimensions of the wound without actually filling the site with new tissue. This phenomenon is called contracture. This type of wound closure can be counteracted by physical splinting of the wound with retaining rings, biopolymer (e.g., collagen) plugs, or other mechanical devices. This rodent model for wound healing with CorMatrix has not been explored and could possibly be comparable to the rabbit model described in this protocol. This should be explored.

REFINEMENT: Since the last IACUC approval, have any study refinements been implemented to reduce the degree of pain or distress experienced by study animals, or have animals of lower phylogenetic status or sentience been identified as potential study/training models in this protocol?

A thin layer of EMLA cream (2.5% Lidocaine/2.5% Prilocaine) can be applied topically with a sterile cotton tip applicator to ear wounds post operatively to help alleviate pain.

REDUCTION: Since the last IACUC approval, have any methods been identified to reduce the number of live animals used in this protocol?

No

9. PUBLICATIONS / PRESENTATIONS: (List any scientific publications and/or presentations that have resulted from this protocol. Include pending/scheduled publications or presentations).

None at this time

10. Were the protocol objectives met, and how will the outcome or training benefit the DoD/USAF?

No, the original objectives of this protocol (pilot study) were to compare healing over time, scar formation and establish a library of histopathology slides for wounds with and without CorMatrix®. Unfortunately, technique perfecting was more challenging than anticipated and the protocol was amended (with IACUC approval) to use 5 animals for technique perfecting instead of only 1. While this increase in animal number did not exceed the original number of 10 animals approved for this protocol, it did not leave enough remaining animals available to conduct the actual pilot study. A new protocol will be submitted to continue the work.

Valuable information was gained from the technique perfecting work that will help investigators move forward with this wound healing model.

11. PROTOCOL OUTCOME SUMMARY: (Please provide, in "ABSTRACT" format, a summary of the protocol objectives, materials and methods, results - include tables/figures, and conclusions/applications.)

Objective: Technique perfecting for wound healing model in the rabbit

Results: Using the five rabbits for technique perfecting proved very beneficial for the investigator. Multiple lessons were learned.

1. The CorMatrix covered the ear wounds best when they were made slightly larger than the wound. Originally, the ear wounds and CorMatrix patches were both cut with a 6mm biopsy punch. This made stitching the patch in place on all sides difficult. The patch was flush with the wound on one side and not necessarily on the opposite side. When the wound was made with a 5mm biopsy punch and the patch with the 6mm punch, the patch easily covered the entire wound on all sides when secured in place.
2. Four wounds per ear was optimum. More than four wound made it difficult (if not impossible) to avoid lacerating a large blood vessel or creating a wound proximal to the ear orifice (an area difficult to place a wound dressing). It also made it difficult to maintain consistent spacing between each wound in an area of the ear that could have a dressing applied.
3. The original spacing between the wounds was 30 mm. That is much larger than necessary or desired. The 30mm distance made applying a dressing to cover all four wounds challenging. Making the wounds 10-15mm apart kept the wounds far enough apart to easily suture the treatment patches in place and not cross contaminate the controls. It also made it much easier to apply a dressing to the ears.
4. Each stitch required a minimum of four throws to tack the CorMatrix in place on four quadrants. The PDS suture has memory and easily unties when not adequately secured. In the first animals, only two throws per stitch were performed and the patches did not remain in place. Upon examination of corresponding histology slides there was no evidence of the CorMatrix patches. When four throws were used to secure the patches, the patch was verified by the histology results.
5. The 5-0 PDS worked much better than 5-0 Vicryl for suturing the patches in place. The Vicryl did not move through the tissue and cartilage smoothly (it was very rough).
6. The Tegaderm dressing worked best when only a small piece was (cut to size) used on the ear. However it would only last one day at best, before coming off the ear. Bacitracin ointment applied twice a day worked well to keep the wound beds moist. However, the bacitracin ointment applied to the wounds with a Vaseline/petroleum gauze placed over it appeared to be the best ear dressing (for maintain wound moisture) and would remain in place for at least a day (possibly two days) before changing.
7. EMLA cream could also be applied post operatively before the dressing is applied for postoperative analgesia.
8. The wounds healed very quickly and were completely healed before 28 days post wounding.



(PI / TC Signature)

15 May 2014
(Date)